



Soft Contact Lens Care Guide

BAUSCH + LOMB

See better. Live better.

Est. 1853



Congratulations.


**You have been fitted with your new
Bausch + Lomb soft contact lenses.**

At Bausch + Lomb our mission is to protect and enhance the precious gift of sight for every individual. We are dedicated to providing a breadth of outstanding products that will dramatically improve the vision and lives of billions of people. As a healthcare company, we hold ourselves to the highest standards of ethics, quality and safety.

Working with our partners and practitioners, we pioneer new technologies and discoveries that continue to bring new vision to the eyes of the world.

To get the best out of your lenses, you need to look after them. This handbook will show you exactly how to do that.

If you need further information, please discuss any concerns you have with your optician.



**I
N
D
E
X**

+ How do I use my lenses correctly?

- Checking your lenses
- Applying your lenses
- Removing your lenses

+ How do I take care of my lenses?

- When applying your lenses
- When removing your lenses

+ Make sure that you DO

+ Make sure that you DON'T

+ Frequently Asked Questions

BAUSCH + LOMB

See better. Live better.


How do I use my lenses correctly?

Applying and removing lenses

Here are a few things to remember that will soon become second nature. It's mostly a matter of routine and a bit of hand/eye co-ordination. Do consult your optician if you are having any difficulty with application or removal.

First, check your lens:

- + Wash your hands with a mild soap and dry them with a lint-free towel.
- + Avoid mix-ups by always putting the same lens in first.
- + Place your lens on the tip of the index finger and check it is not damaged.
- + Check the edges to make sure that the lens is the right side out.
- + If the corners turn out, it's inside out. So turn it the right way and place it back on the tip of your index finger.



Applying your lens:

- + Look straight ahead and pull down your lower eyelid with the middle finger of the same hand.
- + Gently pull the top lid up with the middle finger of the other hand and look up.
- + Place the lens on your eye.
- + Look up, down, left and right to expel any air bubbles.
- + Remove your middle finger, release your eyelid and look down.
- + Then blink a couple of times to make sure it's in place.
- + Now, repeat with the other eye.

Removing your lens:

- + Wash your hands with a mild soap and dry them with a lint-free towel.
- + Now look up and pull your lower lid down with the middle finger.
- + Place your index finger on the lower edge of your lens and slide the lens down to the white of your eye.
- + Lightly squeeze the lens between your index finger and thumb and remove gently.
- + Now, repeat with the other eye.

BAUSCH + LOMB

See better. Live better.

How do I take care of my lenses?

The key to maintaining healthy eyes is to care for your lenses with the best cleaning and care regime. Bausch + Lomb is committed to raising awareness and to providing education on healthy lens care, so we can help you achieve the best results from your lenses.

As part of a healthy lens care routine, follow these simple steps:*

When applying your lenses



Wash your hands with mild soap and dry them with a lint-free towel.



Remove your contact lens from the case / blister pack. Check the lens is the right way out and undamaged.



Place lens on eye.



For reusable lenses, make sure you empty, rinse and wipe your lens case using fresh contact lens solution and a clean tissue.*



Place the case and caps face down on a clean tissue.*

* Please note that this care routine is not necessary for daily disposable lens users.

**Always refer to the instructions for your prescribed solution.



When removing your lenses



Wash your hands with mild soap and dry them with a lint-free towel.



Fill your clean lens case just below the rim with contact lens solution.



Remove your contact lens. Place the lens in the palm of your hand. Apply 3 drops of lens solution to each side of the lens and gently rub for 20 seconds on both sides.



Thoroughly rinse each side of the lens for 5 seconds with fresh lens solution.



Place cleaned lens in the lens case and fill with fresh lens solution. Soak for at least 4-6 hours.**

Ask yourself these three questions, each time you wear your lenses:

- 1 Do my eyes feel good with my lenses? **No discomfort**
- 2 Do my eyes look good? **No redness**
- 3 Do I see well? **No unusual blurring with either eye**

If the answer to any of these questions is no, stop wearing your lenses and consult your contact lens practitioner immediately, who will advise you on what to do **next**.

BAUSCH + LOMB
See better. Live better.

A clean lens means a comfortable lens

Make sure that you DO:

- ✔ Wash your hands and dry them thoroughly before touching your lenses or your eyes
- ✔ Follow the lens care advice as given by your contact lens practitioner and the instructions with the solution
- ✔ Always use fresh solution - never re-use solution
- ✔ Clean your lens case daily with your lens care solution; wipe out your lens case with a clean tissue and leave upside down on a fresh tissue to air dry
- ✔ Replace your lens case as recommended by your contact lens practitioner - usually monthly
- ✔ Follow the directions of your contact lens practitioner for proper use of contact lenses and lens care products

* This is not an exhaustive list that covers all aspects for contact lens usage.





Make sure that you DON'T:

- ✘ Use tap water, or any other water, on your lenses or lens case
- ✘ Wear your lenses if your eyes are red or irritated
- ✘ Sleep or nap in your lenses unless indicated for overnight wear
- ✘ Touch the bottle tip to any surface or your eye which could cause contamination to the eye or solution
- ✘ Apply drops or eye medication without first consulting your contact lens practitioner, doctor or pharmacist
- ✘ Top up the solution in your lens case or decant your solution into another container
- ✘ Switch your solution without the advice of your contact lens practitioner
- ✘ Re-use daily disposables. They are for single-use only.
- ✘ Wear contact lenses while showering, bathing or swimming

BAUSCH + LOMB

See better. Live better.

FAQs

What happens if I keep wearing my lenses after I should have thrown them out?

It is important to replace your lenses as advised by your contact lens practitioner. Even over just a short time, microscopic deposits and bacteria are formed on the lens surface, which means they become less comfortable and carry a higher risk of infection.


What happens if I sleep or nap with my lenses in?

Contact lenses should be removed before sleeping.* Sleeping in lenses places you at a higher risk of suffering serious eye infections.

* Unless you are using specialist lenses designed for overnight wear, upon the recommendation of your contact lens practitioner.

Why shouldn't I use tap water to clean and store my contact lenses?

Tap water carries micro-organisms that can cause serious eye infections. So keep it away from your contact lenses and lens case. Ensure your contact lenses and lens case are cleaned as recommended by your contact lens practitioner.



Give your contact lens practitioner a look-in

Now that your Contact Lens Practitioner has prescribed the correct Bausch + Lomb contact lenses for your needs, they will advise you on renewing your prescription and ordering replacements. They will also advise you on the regularity of your eye examinations and contact lens check-ups.

My Contact Lens

Right lens

.....

Left lens

.....

Phone number of my contact lens practitioner

.....

My next check-ups

1.


3.

2.

4.

BAUSCH + LOMB

See better. Live better.



www.bausch.co.uk