

Bio true



ONEday lenses

for PRESBYOPIA



Give your patients stable vision across all distances with Biotrue® ONEday lenses for Presbyopia

Parameters

Material	nesofilcon A
Water content	78%
Oxygen	42 DK/t at centre -3.00 allowing 93% oxygen flux
Powers	-9.00 to +6.00D (0.25 steps) Spectacle Add Powers: Low: up to +1.50D High: from +1.75D
Design	3-Zone Progressive Design™ 
Base curve	8.6mm
Diameter	14.2mm
Centre thickness	0.1mm for -3.00
Wear	Daily
UV protection*	Yes 
Packs	30 pack



Please contact your Bausch & Lomb Territory Manager for more information

order online: www.bauschonline.co.uk

UK: 0845 602 2350

UK Fax: 0845 602 2351

ROI: 1800 409077

ROI Fax: 1800 409083

BAUSCH + LOMB

See better. Live better.

*WARNING: UV-absorbing contact lenses are NOT substitutes for protective UV-absorbing eyewear, such as UV-absorbing goggles or sunglasses, because they do not completely cover the eye and surrounding area. The effectiveness of wearing UV-absorbing contact lenses in preventing or reducing the incidence of ocular disorders associated with exposure to UV light has not been established at this time. You should continue to use UV-absorbing eyewear as directed. NOTE: Long-term exposure to UV radiation is one of the risk factors associated with cataracts. Exposure is based on a number of factors such as environmental conditions (altitude, geography, cloud cover) and personal factors (extent and nature of outdoor activities). UV-blocking contact lenses help provide protection against harmful UV radiation. However, clinical studies have not been done to demonstrate that wearing UV-blocking contact lenses reduces the risk of developing cataracts or other eye disorders. Biotrue and inspired by the biology of your eyes are trademarks of Bausch & Lomb Incorporated or its affiliates. All other product/brand names are trademarks of their respective owners. © 2015 Bausch & Lomb Incorporated. US/BPR/15/0002

How to fit Biotrue® ONEday for Presbyopia

Select initial lenses

- Update spectacle refraction and Add power
- Determine ocular dominance for distance vision
- Select contact lens distance prescription based upon spherical equivalent from spectacle Rx, adjusted for vertex distance if necessary
- Choose trial contact lenses based upon the above calculation and select add power

ADD Selection

	SPECTACLE ADD	BOTH EYES
UP TO	<1.50	Low Add
OVER	>1.50	High Add

Suggested Patient Criteria

- Good motivation and realistic expectations
- Refractive astigmatism no greater than -1.00D

Evaluate initial lenses

- Allow trial lenses to equilibrate for at least 10 minutes before assessing fit and vision
- Evaluate distance and near vision binocularly in normal room illumination
- If vision at distance and near are satisfactory, dispense lenses and schedule follow-up exam within 1-2 weeks

To refine near vision

If patient is wearing two Low Add lenses:

	DOMINANT EYE	NON-DOMINANT EYE
INITIAL LENS	Low Add	Low Add
REFINEMENT 1	Low Add	High Add

Refinement 2: If vision is still unsatisfactory, make small changes by adding +0.25D at a time to non-dominant eye (wearing High Add lens) using handheld lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.

If patient is wearing two High Add lenses:

	DOMINANT EYE	NON-DOMINANT EYE
INITIAL LENS	High Add	High Add
REFINEMENT 1	High Add	Add +0.25D to the non-dominant eye

Refinement 2: If vision is still unsatisfactory, make small changes by adding +0.25D at a time to non-dominant eye using handheld lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.

To refine distance vision

If patient is wearing two Low Add lenses:

	DOMINANT EYE	NON-DOMINANT EYE
INITIAL LENS	Low Add	Low Add
REFINEMENT 1	Fit Biotrue® ONEday Sphere	Low Add

Refinement 2: If vision is still unsatisfactory, make small changes by adding -0.25D at a time to dominant eye (wearing Biotrue® ONEday spherical lens) using handheld lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.

If patient is wearing two High Add lenses:

	DOMINANT EYE	NON-DOMINANT EYE
INITIAL LENS	High Add	High Add
REFINEMENT 1	Low Add	High Add

Refinement 2: If vision is still unsatisfactory, make small changes by adding -0.25D at a time to dominant eye (wearing Low Add lens) using handheld lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.

Biotrue® ONEday for Presbyopia has more consistent 'Add' power

across the available power range than Dailies AquaComfort Plus Multifocal Low Add and Proclear 1 Day Multifocal.

