



INSTRUCTIONS FOR USE

This leaflet is intended to inform users of Bausch + Lomb daily disposable hydrophilic soft contact lenses how to wear and care for their lenses.

It applies to:

PRODUCT NAME	MATERIAL	WATER CONTENT	PACKAGING SOLUTION
SofLens® daily disposable (hilafilcon B) Visibility Tinted Contact Lenses	Hilafilcon B	59% by weight	Borate buffered saline with Poloxamine
SofLens® daily disposable Toric (hilafilcon B) Visibility Tinted Contact Lenses			
Naturelle™ daily disposable (hilafilcon B) Cosmetically Tinted Contact Lenses			
Biotrue® ONEday (nesofilcon A) Visibility Tinted Contact Lenses	Nesofilcon A	78% by weight	
Biotrue® ONEday For Presbyopia (nesofilcon A) Visibility Tinted Contact Lenses			
Biotrue® ONEday For Astigmatism (nesofilcon A) Visibility Tinted Contact Lenses			

Indication

The Bausch + Lomb (hilafilcon B and nesofilcon A) Contact Lens is indicated for the daily wear correction of refractive ametropia (myopia, hyperopia, astigmatism), presbyopia and reduction of spherical aberration in aphakic and/or non-aphakic persons with non-diseased eyes that exhibit refractive astigmatism up to 2.00 diopters or less (SofLens® daily disposable Toric and Biotrue® ONEday For Astigmatism are 5.00 Diopters or less), that does not interfere with visual acuity. The lens is to be prescribed for single-use disposable wear in spherical powers ranging from +20.00D to -20.00D, and is to be discarded after each removal.

Description and wear application

These lenses are intended for daily disposal wear and should be discarded after each wearing period. Reuse without cleaning and disinfecting the lenses could lead to infection and/or serious eye injury. Always follow your eye care practitioner's recommended wearing schedule.




These lenses are supplied sterile and are intended by the manufacturer for a maximum wear period of up to one day. They have been sterilized by steam.

If the packaging is damaged or opened DO NOT USE.

Insertion and removal instructions

Before handling lenses, always wash and dry hands thoroughly avoiding perfumed or medicated soaps. Pour the lens into the palm of your hand. Make sure the lens is the right way out.



Inspect your lens before inserting and do not insert if scratched or damaged.

To remove your lens, make sure your hands are clean and dry, use your finger to slide the lens down to the white part of the eye and remove it. Use a rewetting drop if the lens is difficult to remove. Always discard lens following removal from the eye.

Caring for your lenses

- Always wash and rinse hands before handling lenses
- Do not allow cosmetics, lotions, soaps, creams or perfumed substances to get in the eyes or on lenses
- Never use sharp objects or finger nails to remove lenses from their containers or your eyes
- Always discard lenses after removal

Warnings and precautions

Failure to follow your eye care practitioner's recommended wear schedule and cleaning regime can give rise to greater risk of serious eye infections.

Do not wear daily disposable lenses for longer than a day or while sleeping.

Do not wear contact lenses if you are suffering from:

- Allergy, inflammation, infection or redness in or around the eye
- Dry eyes (inadequate tear fluid)
- Poor health affecting the eyes, e.g. colds and flu
- Any systemic illness affecting the eyes

Contact your eye care practitioner before wearing contact lenses when:

- Practicing water sports
- In excessively dry or dusty environments
- Using eye medication

If you experience any of the following problems, remove your lenses immediately:

- Eyes stinging, burning, itching, irritation or other eye pain
- Reduced comfort compared with previous wearing
- Unusual secretions or excessive watering from the eye
- Redness of eyes
- Severe or persistent dry eyes
- Reduced or blurred vision
- Rainbows or halos around lights
- Sensitivity to light

If any of these symptoms continue after removal and insertion of a new lens, then remove the lens and contact your eye care practitioner. These symptoms can lead to more serious conditions if ignored.

Additional warnings and precautions for Biotrue ONEday (nesofilcon A) Contact Lens

Warning: UV absorbing contact lenses are NOT substitutes for protective UV absorbing eyewear such as UV absorbing goggles or sunglasses because they do not completely cover the eye and surrounding area. You should continue to use UV absorbing eyewear as directed.

Note: Long term exposure to UV radiation is one of the risk factors associated with cataracts. Exposure is based on a number of factors such as environmental conditions (altitude, geography, cloud cover) and personal factors (extent and nature of outdoor activities). UV blocking contact lenses help provide protection against harmful UV radiation.

Note: The effectiveness of wearing UV-absorbing contact lenses in preventing or reducing the incidence of ocular disorders associated with exposure to UV-light has not been established at this time. However, clinical studies have not been done to demonstrate that wearing UV blocking contact lenses reduce the risk of developing cataracts or other eye disorders. Consult your Eye Care Professional for more information.

Reference guide for symbols, labels and cartons



Fee paid for waste management per 94/62/EC



Bausch & Lomb Incorporated
1400 North Goodman Street
Rochester NY, 14609, USA



Bausch & Lomb Incorporated
Cork Road Industrial Estate
Waterford, X91 V383, Ireland

Date of Issue: 2018-06



0050 Indicates the CE conformity marking and the notified body number



Do Not Reuse



Lower limit of temperature
1°C / 34°F



See Instruction for Use



Sterile using steam or dry heat



For US only: Caution Federal law restricts this device to sale by or on the order of a licensed practitioner

®/™ are trademarks of Bausch & Lomb Incorporated or its affiliates.

© Bausch & Lomb Incorporated.

1802322
P9EU05