



RANGE  
EXPANDED  
NOW 3280  
PARAMETERS

## SIMPLE STEPS FOR FITTING

- STEP 1:** Update spectacle refraction and Add power  
Determine the dominant eye using a sensory dominance test with a handheld +1.50 trial lens
- STEP 2:** Select toric lens power (adjust for vertex distance if necessary)
- STEP 3:** Select **Add power** following guidance below:

## EVALUATE THE LENS FOR SUCCESS

- Allow trial lenses to equilibrate for at least 15 minutes before assessing fit and vision
- Confirm axis orientation
- Evaluate distance and near vision binocularly in normal room illumination
- If vision at distance and near are satisfactory, dispense lenses and schedule follow-up exam within 1-2 weeks

## RECOMMENDATIONS

Set realistic patient expectations in a positive manner

Use Bausch + Lomb fitting guide to achieve the best results

Use most up to date spectacle refraction to select initial lens

**DO NOT** use the current prescription for other multifocal lenses

## REFINE IF NEEDED: FOLLOW GUIDANCE BELOW

NEAR VISION			DISTANCE VISION		
	DOMINANT EYE	NON-DOMINANT EYE		DOMINANT EYE	NON-DOMINANT EYE
TWO LOW ADDS	INITIAL LENS	LOW ADD	INITIAL LENS	LOW ADD	LOW ADD
	REFINEMENT 1	LOW ADD	REFINEMENT 1	Bausch + Lomb ULTRA® for Astigmatism	LOW ADD
<b>Refinement 2:</b> If vision is still unsatisfactory, make small changes by adding +0.25D at a time to non-dominant eye (wearing High Add lens) using hand-held lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.			<b>Refinement 2:</b> If vision is still unsatisfactory, make small changes by adding -0.25D at a time to dominant eye (Bausch + Lomb ULTRA® for Astigmatism) using hand-held lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.		
TWO HIGH ADDS	INITIAL LENS	HIGH ADD	INITIAL LENS	HIGH ADD	HIGH ADD
	REFINEMENT 1	HIGH ADD	REFINEMENT 1	LOW ADD	HIGH ADD
<b>Refinement 2:</b> If vision is still unsatisfactory, make small changes by adding +0.25D at a time to non-dominant eye using hand-held lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.			<b>Refinement 2:</b> If vision is still unsatisfactory, make small changes by adding -0.25D at a time to dominant eye (wearing Low Add lens) using hand-held lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.		

## PARAMETERS

MATERIAL:	sarnfilcon A
LENS MATERIAL TECHNOLOGY:	MoistureSeal® technology
WATER CONTENT:	46%
LENS DESIGN TECHNOLOGY:	3-Zone Progressive Design™, OpticAlign™ Design
BASE CURVE:	8.6 mm
DIAMETER:	14.5 mm
SPHERICAL POWERS:	+4.00D to -6.00D in 0.25D steps including plano
CYLINDER POWERS:	-0.75D, -1.25D, -1.75D and -2.25D
AXES:	10°, 20°, 70°, 80°, 90°, 100°, 110°, 160°, 170°, 180°
ADD POWERS:	Low: Up to +1.50D   High: From +1.75D
VISIBILITY TINT:	Light blue
WEAR SCHEDULE:	Daily wear, up to 7 days extended wear
REPLACEMENT SCHEDULE	Monthly

## ADD SELECTION

SPECTACLE ADD	BOTHEYES
Up to +1.50D	LOW ADD
From +1.75D	HIGH ADD



## SIMPLE STEPS FOR FITTING

- STEP 1:** Update spectacle refraction and Add power  
Determine the dominant eye using a sensory dominance test with a handheld +1.50 trial lens
- STEP 2:** Select contact lens distance prescription based upon spherical equivalent from spectacle Rx and following (adjusted for vertex distance if necessary)
- STEP 3:** Select **Add power** following guidance below:

## EVALUATE THE LENS FOR SUCCESS

- Allow trial lenses to equilibrate for at least 15 minutes before assessing fit and vision
- Evaluate distance and near vision binocularly in normal room illumination
- If vision at distance and near are satisfactory, dispense lenses and schedule follow-up exam within 1-2 weeks

## RECOMMENDATIONS

Set realistic patient expectations in a positive manner

Use Bausch + Lomb fitting guide to achieve the best results

Use most up to date spectacle refraction to select initial lens

**DO NOT** use the current prescription for other multifocal lenses

## REFINE IF NEEDED: FOLLOW GUIDANCE BELOW

NEAR VISION			DISTANCE VISION		
	DOMINANT EYE	NON-DOMINANT EYE		DOMINANT EYE	NON-DOMINANT EYE
IF PATIENT IS WEARING:	TWO LOW ADDS		INITIAL LENS	LOW ADD	LOW ADD
			REFINEMENT 1	LOW ADD	HIGH ADD
<p><b>Refinement 2:</b> If vision is still unsatisfactory, make small changes by adding +0.25D at a time to non-dominant eye (wearing High Add lens) using handheld lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.</p>			<p><b>Refinement 2:</b> If vision is still unsatisfactory, make small changes by adding -0.25D at a time to dominant eye (Bausch + Lomb ULTRA® spherical lens) using hand-held lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.</p>		
TWO HIGH ADDS	TWO HIGH ADDS		INITIAL LENS	HIGH ADD	HIGH ADD
			REFINEMENT 1	HIGH ADD	ADD +0.25D TO THE NON-DOMINANT EYE
<p><b>Refinement 2:</b> If vision is still unsatisfactory, make small changes by adding +0.25D at a time to non-dominant eye using hand-held lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.</p>			<p><b>Refinement 2:</b> If vision is still unsatisfactory, make small changes by adding -0.25D at a time to dominant eye (wearing Low Add lens) using hand-held lenses, and continue evaluating vision binocularly in normal room illumination. Adjust contact lens power when vision is satisfactory.</p>		

## PARAMETERS

MATERIAL:	silfilcon A
LENS MATERIAL TECHNOLOGY:	MoistureSeal® Technology
WATER CONTENT:	46%
LENS DESIGN TECHNOLOGY:	3-Zone Progressive Design™ with centre-near aspheric optics
BASE CURVE:	8.5 mm
DIAMETER:	14.2 mm
SPHERICAL POWERS:	+6.00D to -10.00D (in 0.25D steps including plano)
ADD POWERS:	Low: up to +1.50D   High: From +1.75D
VISIBILITY TINT:	Light blue
WEAR SCHEDULE:	Daily wear, up to 7 days extended wear
REPLACEMENT SCHEDULE	Monthly

## ADD SELECTION

SPECTACLE ADD	BOTH EYES
Up to +1.50D	LOW ADD
From +1.75D	HIGH ADD